











Eye Gaze Setup Guide

Considerations and tips for setting up Eye Gaze with your learners to get the best results.

| | |
|---|---|
| <p>1 Positioning Equipment</p>  <ul style="list-style-type: none">✓ Adjust screen height with appropriate stand.✓ Adjust seating if possible.✓ Use height adjustable table if required. | <p>2 Positioning Learner</p>  <ul style="list-style-type: none">✓ Adjust eye tracker/screen angle to suit learner.✓ Learner should be 50-75cm away from screen.✓ Make sure eye tracker is level with or just below learner's eyes. |
| <p>3 Reclined Position</p>  <ul style="list-style-type: none">✓ If learner is usually in a reclined position, mount screen above using specialist floor stand. | <p>4 Tilted Position</p>  <ul style="list-style-type: none">✓ If learner's head is normally tilted, you can tilt the screen/eye tracker to match. |
| <p>5 Lighting</p>  <ul style="list-style-type: none">✓ Both the learner and the eye tracker will struggle to see if conditions are too bright e.g. strong overhead lights or bright sunlight. | <p>6 Glasses</p>  <ul style="list-style-type: none">✓ Eye trackers work with most glasses. Some specialist lenses e.g. varifocals might not work so well. |
| <p>7 Track both or just 1 eye</p>  <ul style="list-style-type: none">✓ If learner has significant visual problem in one eye, set eye tracking for the other eye only. <p>Activate this setting from within the eye tracker's setup software.</p> | <p>8 Eye Track status</p>  <ul style="list-style-type: none">✓ Always check learner's eyes are in a central position in track status window before you start. Depending on the eye tracker, the track status window can be accessed from within the eye tracker's setup software, from a keyboard shortcut or as a 'red amber green' on screen visual guide. |
| <p>9 Calibration</p>  <ul style="list-style-type: none">✓ Calibrate the learner's eyes within the eye tracker setup software and save their User Profile. You do not have to re-calibrate every time they use the eye tracker. Just check their User Profile is selected. Individual points can be improved if required. | <p>10 Windows Control</p>  <ul style="list-style-type: none">✓ If software being used is not 'eye-gaze enabled' make sure Windows Control is on within the eye tracker setup software and that the correct mouse setting is applied (usually 'single left click'). |

For more information on Eye Gaze Technology go to: www.callscotland.org.uk/information/eye-gaze/